

Flanges

Evactron Adaptor Flanges –SEM Port to KF40

Part Number

Description

OD: Outside diameter of flange

ORID: O Ring Inside Diameter O-rings not included

H: mounting holes #

Dimensions are in millimeters

FEI

230090-01.....FLANGE.....FEI GIS ASYM HEX,
230141-01.....FLANGE.....FEI OD-100 ORID-75, 6H symmetric
230153-01.....FLANGE.....FEI OD 85 ORID 60, 3H
230154-01.....FLANGE.....FEI OD-100 ORID-70, 3H
230235-01.....FLANGE.....FEI TECHNAI PLUG
230335-01.....FLANGE.....FEI 235 OD 88, ORID 60, 3H
230350-01.....FLANGE.....FEI LARGE OD 120, ORID 60,
230351-01.....FLANGE.....FEI Stata Y port I PIECE OD 59 ORID 38
230359-01.....FLANGE.....FEI OD 90 ORID 70 3H

HITACHI

230568-01.....FLANGE.....HITACHI OD-69, ORID-47.5, 4H SYM
230320-01.....FLANGE.....HITACHI S 5000 PLUG, 4H
230333-01.....FLANGE.....HITACHI S5200 PLUG, 4H

ISI

230588-01.....FLANGE.....ISI RECT 330 X180 4H

JEOL

230001-01.....FLANGE.....JEOL 6700f, OD 64 ORID 40, 4H
230002-01.....FLANGE.....JEOL 99 OD, 75 ORID, 4H
230311-01.....FLANGE.....JEOL 6500, OD 88, ORID, 45, 4 ASYM
230322-01.....FLANGE.....JEOL 6320 Arched 98HX75W 4H
230323-01.....FLANGE.....JEOL 5900, OD 94, ORID 65 4H
230324-01.....FLANGE.....JEOL 6360, OD 99, ORID 75 4H
230664-01.....FLANGE.....JEOL OM port OD-95, ORID-76, 3H SYM
230672-01.....FLANGE.....JEOL OD 160, ORID 110
230893-01.....FLANGE.....JEOL JSM-35,WDS port

LIECA

230012-01.....FLANGE.....LEICA 360, 64SQ 4H

ZEISS (LEO)

230019-01.....FLANGE.....ZEISS OD 78, ORID 53, 4H
230034-01.....FLANGE.....ZEISS LEFT 96mm, 72 ORID 4H sym
230356-01.....FLANGE.....ZEISS LF 100, ORID 60, 3H
230357-01.....FLANGE.....ZEISS OD 89, ORID 60, 4H ASYM
230358-01.....FLANGE.....ZEISS EDS, RECT 97 x 81
230364-01.....FLANGE.....ZEISS SQ 85MM ORID 70
230366-01.....FLANGE.....ZEISS WDS 190x150

XEI SCIENTIFIC, INC.
1755 EAST BAYSHORE
SUITE 17
REDWOOD CITY,
CA 94063
MAIN 650.369.0133
800.500.0133
FAX 650.363.1659
HTTP://EVACTION.COM
INFO@EVACTION.COM
SALES@EVACTION.COM
1.415.566.5774



XEI

SCIENTIFIC, INC.

